



HSCo
WEIGHING SOLUTIONS



DIGITAL DISPLAY MAGNETIC STIRRERS

PRO SERIES

Precision Mixing. Reliable Results.



POWERFUL MAGNETIC STIRRING

Strong magnetic drive for consistent and efficient mixing.



DIGITAL DISPLAY

Clear screen for accurate speed and performance monitoring.



PRECISE SPEED CONTROL

Adjustable speed control knob for smooth and precise operation.



DURABLE & CHEMICAL RESISTANT

Built with high-quality materials for long-lasting performance.



SAFE & USER FRIENDLY

Stable design with non-slip base and easy-to-use interface.



PRODUCT OVERVIEW

The HSCo Digital Display Magnetic Stirrer – Pro Series is engineered for precise and efficient mixing in laboratories. It features a strong magnetic drive, digital speed display, and smooth speed control for reliable performance in various applications. Its durable design ensures safety, stability, and long-lasting use.

IDEAL FOR



Chemical Labs



Research & Development



Pharmaceutical Labs



Educational Institutions



Quality Control Labs

TECHNICAL SPECIFICATIONS

Model	: Digital Display Magnetic Stirrer – Pro Series
Stirring Capacity (H ₂ O)	: Up to 5 Liters
Speed Range	: 200 – 1500 RPM
Speed Display	: Digital
Speed Control	: Knob Control
Stirring Bar Length	: Up to 50 mm
Top Plate Material	: ABS Plastic (Chemical Resistant)
Motor Type	: Brushless DC Motor
Power Supply	: 220 – 240V, 50/60Hz
Power Consumption	: 20W
Operating Temperature	: 5°C to 40°C
Permissible Relative Humidity	: 80%
Dimensions (L x W x H)	: 210 x 150 x 70 mm
Weight	: Approx. 0.85 kg

WHY CHOOSE HSCo MAGNETIC STIRRERS?



Consistent Performance



High Precision & Accuracy



Durable Construction



Easy to Use & Maintain



Value for Money

CONTACT US



Mobile no.:
7045922250



Mail:
ask@hindustanscale.com



Website :
www.hindustanscale.com