



**HSCo**  
WEIGHING SOLUTIONS

# BMI 250P Heavy Duty – Weighing Gym Height Scale



High  
Precision  
Sensors



User  
Friendly  
Platform



Compact  
Design



7045922250



ask@hindustanscale.com



www.hindustanscale.com



# BMI 250P Heavy Duty - Weighing Gym Height Scale



## PRODUCT OVERVIEW

The BMI 250P Heavy Duty Weighing Gym Height Scale is a professional-grade health monitoring device designed for gyms, hospitals, clinics, and diagnostic facilities, providing accurate weight, height, and BMI measurements with robust and durable construction.



## CORE FEATURES



### High Precision Weighing

Advanced load-cell technology ensures accurate weight measurement



### Integrated Height Measurement

Built-in height measuring system for quick height assessment




### Automatic BMI Calculation

Calculates BMI based on height for instant health assessment



### Heavy Duty Construction

Strong and durable body ideal for intensive usage in gyms, hospitals, and fitness centers



### Digital Display Interface

Clear digital display provides easy-to-read results for weight, height, and BMI metrics



### Professional Healthcare Use

Ideal for clinical health monitoring. & routine physical examinations



7045922250



ask@hindustanscale.com



www.hindustanscale.com

# TECHNICAL SPECIFICATIONS



**1** Ultrasonic probe height measurement

**2** Folding height Rod

**3** High-definition LCD digital display screen

**4** Printing function:  
With thermal printer

**5** Column

**6** Weighing platform

Voice reading function  
(switchable)

Switch on/off mode:  
Gravity/button on, Button off

Weight measurement:  
3-200 KG

Weight Accuracy:  
50 G

Height Measurement:  
50 - 200 CM

Height Accuracy:  
0.5 CM

Weight unit :  
Kg/LB

Display screen :  
70\*80mm VA screen

With RS232 port:  
Data output can be realized

Function  
Weight, height, BMI, voice function,  
data record, review function